

NELFILCON A DAILY DISPOSABLE SOFT CONTACT LENSES



INSTRUCTIONS FOR USE - THIS LEAFLET CONTAINS IMPORTANT PRODUCT USE AND SAFETY INFORMATION. PLEASE READ CAREFULLY AND RETAIN FOR FUTURE REFERENCE. CONTACT LENSES SHOULD ALWAYS BE FITTED BY AN EYE CARE PROFESSIONAL. IT IS ESSENTIAL THAT YOU FOLLOW YOUR EYE CARE PROFESSIONAL'S DIRECTIONS AND ALL LABELLING INSTRUCTIONS FOR PROPER USE OF YOUR CONTACT LENSES.



CAUTION: Federal (USA) law restricts this device to the sale by or on the order of a licensed eye care professional

PRODUCT NAMES

DAILIES® AquaComfort Plus®
DAILIES® AquaComfort Plus® Toric
DAILIES® AquaComfort Plus® Multifocal
DAILIES® FRESHLOOK® ILLUMINATE
FOCUS® DAILIES®
FOCUS® DAILIES® Progressives
FOCUS® DAILIES® Toric
FRESHLOOK® One-Day

MATERIAL AND PACKAGING

DAILIES® AquaComfort Plus® soft contact lenses / **DAILIES® AquaComfort Plus® Multifocal** soft contact lenses / **DAILIES® AquaComfort Plus® Toric** soft contact lenses contain 31% nelfilcon A and 69% water.

The lenses are tinted light blue to assist handling; the tint contains Cu-phthalocyanine.

DAILIES® AquaComfort Plus® / DAILIES® AquaComfort Plus® Multifocal / DAILIES® AquaComfort Plus® Toric are packaged in strips of 5 individual blisters containing isotonic, phosphate-acetate buffered saline. The saline contains HPMC and PEG400. The saline may contain up to 0.05 % Poloxamer.

DAILIES® FRESHLOOK® ILLUMINATE soft contact lenses contain 31% nelfilcon A and 69% water. The lenses are printed with an intermittent coating containing a combination of the following pigments: iron oxides, carbazole violet, PCN blue, titanium dioxide.

DAILIES® FRESHLOOK® ILLUMINATE contact lenses are packaged in strips of 5 individual blisters containing isotonic, phosphate-acetate buffered saline. The saline may contain up to 0.02% Poloxamer.

FOCUS® DAILIES® soft contact lenses / **FOCUS® DAILIES® Progressives** soft contact lenses / **FOCUS® DAILIES® Toric** soft contact lenses contain 31% nelfilcon A and 69% water. The lenses are tinted light blue to assist handling; the tint contains Cu-phthalocyanine.

FOCUS® DAILIES® / FOCUS® DAILIES® Progressives / FOCUS® DAILIES® Toric are packaged in strips of 5 individual blisters containing isotonic, phosphate-acetate buffered saline. The saline may contain up to 0.02% Poloxamer.

FRESHLOOK® One-Day soft contact lenses contain 31% nelfilcon A and 69% water. The lenses are printed with an intermittent coating containing a combination of the following pigments: iron oxides, titanium dioxide, Cu-phthalocyanine, chromium oxide, phthalocyanine green and carbazole violet.

FRESHLOOK® One-Day are packaged in strips of 5 individual blisters containing isotonic, phosphate-acetate buffered saline. The saline may contain up to 0.02% Poloxamer.

The lenses are steam sterilised **STERILE** . Use only if printed foil pack is intact.

INTENDED USE

DAILIES® AquaComfort Plus® contact lenses are intended for daily wear (less than 24 hours while awake) for the optical correction of vision in persons with non-diseased eyes with minimal astigmatism that does not interfere with vision. **DAILIES® AquaComfort Plus® Multifocal** contact lenses are intended for daily

wear (less than 24 hours while awake) for the optical correction of vision including presbyopia in persons with non-diseased eyes with minimal astigmatism that does not interfere with vision.

DAILIES® AquaComfort Plus® Toric contact lenses are intended for daily wear (less than 24 hours while awake) for the optical correction of vision in persons with non-diseased eyes who require additional correction for astigmatism.

DAILIES® FRESHLOOK® ILLUMINATE contact lenses are intended for daily wear (less than 24 hours while awake) in one or both eyes for the optical correction of vision in persons with non-diseased eyes with minimal astigmatism that does not interfere with vision.

The print pattern blends with the natural eye colour for a natural looking enhancement.

FOCUS® DAILIES® contact lenses are intended for daily wear (less than 24 hours while awake) for the optical correction of vision in persons with non-diseased eyes with minimal astigmatism that does not interfere with vision.

FOCUS® DAILIES® Progressives contact lenses are intended for daily wear (less than 24 hours while awake) for the optical correction of vision including presbyopia in persons with non-diseased eyes with minimal astigmatism that does not interfere with vision.

FOCUS® DAILIES® Toric contact lenses are intended for daily wear (less than 24 hours while awake) for the optical correction of vision in persons with non-diseased eyes who require additional correction for astigmatism.

FRESHLOOK® One-Day contact lenses are intended for daily wear (less than 24 hours while awake) in one or both eyes for the optical correction of vision in persons with non-diseased eyes with minimal astigmatism that does not interfere with vision.

The print pattern blends with the natural eye colour to create a natural looking change or enhancement.

All nelfilcon A contact lenses are for single use Daily Disposable Wear only.

RECOMMENDED WEAR AND REPLACEMENT SCHEDULE

Nelfilcon A daily disposable contact lenses are intended to be worn once and then discarded at the end of each wearing period. The maximum daily wearing time should be determined by your eye care professional.

Nelfilcon A daily disposable contact lenses have not been developed for use with contact lens cleaners or disinfection systems. Your lenses should be discarded after use. Start each wearing period with a fresh new lens.

DO NOT SLEEP IN YOUR LENSES. Nelfilcon A daily disposable contact lenses have not been developed for wear while sleeping. You should not wear Nelfilcon A daily disposable contact lenses while sleeping.

HANDLING YOUR LENSES

Your eye care professional should provide you with detailed contact lens insertion and removal instructions.

Always wash and dry your hands before handling lenses. Gently shake the blister pack before opening. Peel back the foil lid. Pour the lens into the palm of your hand or if necessary, gently slide the lens out of the container using your forefinger.

Ensure the lens is not inside out and that you have the appropriate lens for that eye. Inspect your lenses prior to insertion, do not insert if damaged.

When removing your lens make sure that your hands are completely dry. Blink fully several times, then while looking up, slide the lens down onto the white part of the eye. Remove the lens by gently pinching it between the thumb and forefinger. If the lens is difficult to remove, do not pinch the eye tissue. Apply a lubricating and rewetting drop and try again after several minutes.

- Never use tweezers, suction cups, sharp objects or your finger nails to remove your lenses from the lens container or your eyes.


WHEN NOT TO WEAR LENSES

Contact lenses should not be worn in conjunction with certain health or environmental conditions. Conditions that may prevent or interfere with safe contact lens wear include:

- Allergy, inflammation, infection or irritation in or around the eye or eyelids.
- Periods of poor health, such as severe colds or influenza.
- The use of some medications, including eye medication.

- Inadequate tear film (dry eyes).
 - Excessively dry or dusty environments that make contact lens wear uncomfortable.
 - Water sport without the use of goggles.
- Consult your eye care professional specifically regarding these or other conditions.

IMPORTANT THINGS TO REMEMBER

- Device is for single use only. Reuse can increase the risk of infection or fatigue-induced device failure.
- **Do not share** your lenses with anyone as this may spread micro-organisms which could result in serious eye health problems.
- Check your eyes daily to ensure they look well, feel comfortable and your vision is clear.
- It is recommended that contact lens wearers visit their eye care professional regularly as directed.
- Do not use lenses after the expiry date. 
- In conditions of low light, wearers of tinted/coloured lenses may experience reduced vision of low contrast objects. Wearers of progressive lenses may experience reduced vision in any light condition.
- Keep a note of the correct lens power for each eye.
- Before lens insertion, check that the lens power on each foil pack is correct for that eye.
- Always carry spare lenses with you.
- Be careful while using soaps, lotions, creams, cosmetics or deodorants as they can cause irritation if they come into contact with your lenses.
- Insert your lenses **before** applying make-up and remove them **before** taking your make-up off.
- Avoid spraying aerosols such as hair spray near your eyes while wearing your lenses as this may cause irritation.
- Consult your eye care professional about wearing lenses during sporting activities including swimming.
- Remove your lenses when exposed to noxious or irritating vapour.
- Never allow contact lenses to come into contact with non-sterile liquids (including tap water and saliva) as microbial contamination can occur, which may lead to permanent eye damage.
- Inform your employer that you wear contact lenses, especially if your job involves using eye protection equipment.
- If a nelficon A daily disposable contact lens is exposed to air while off the eye it will become dry and brittle. It should be discarded and replaced with a new one.

POSSIBLE PROBLEMS

Although soft contact lenses provide many benefits to the wearer, it is possible that problems might occur and may be first noticed as one or more of the following signs:

- Feeling of something in the eye
- Uncomfortable lens
- Eye redness
- Sensitivity to light
- Burning, stinging, itching or watering eyes
- Reduced sharpness of vision
- Rainbows or halos around lights
- Increased eye secretions
- Discomfort/pain
- Severe or persistent dry eyes

These symptoms, if ignored, may lead to more serious complications.

WHAT TO DO IF A PROBLEM OCCURS

If any of the above signs or symptoms occur:

- **IMMEDIATELY REMOVE THE LENS(ES).**
- Inspect your eye in a mirror looking for signs of irregularity.
- If the discomfort or problem stops, insert a fresh, new lens.
- If the above symptoms continue after insertion of a new lens(es) remove the lens(es) **IMMEDIATELY** and promptly contact your eye care professional. A serious condition such as infection, corneal ulcer or iritis may be present. These conditions could progress rapidly and may lead to permanent loss of vision. Less serious reactions such as abrasions, epithelial staining and bacterial conjunctivitis

must be managed and treated early to avoid complications.

- If a lens tears in your eye, it will feel uncomfortable. Do not panic, as it is impossible to lose a contact lens or part of a contact lens behind the eye. Remove the pieces carefully by pinching them as you would do for normal lens removal. If the lens pieces do not seem to remove easily, do not pinch the eye tissue. Rinse with saline. If this does not help, contact your eye care professional for assistance. Your eye care professional can easily find the lens using fluorescein.
- Occasional dryness may be relieved by blinking fully several times or by the use of lubricating and rewetting drops. If dryness persists consult your eye care professional.

ADDITIONAL INFORMATION

The risk of a possible problem can be reduced by recognising the signs and following the advice, "What to do if a Problem Occurs". Contact lens wear increases the risk of eye infections. Smoking further increases the risk of ulcerative keratitis for contact lens users.



MANUFACTURER:
Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, TX 76134-2099, USA



AUTHORIZED EU REPRESENTATIVE:
Alcon Laboratories (UK) Ltd.
Frimley Business Park, Frimley, Camberley
Surrey, GU16 7SR, United Kingdom

Abbreviations for lens specifications:

DIA	Diameter
BC	Base Curve
PWR	Power
CYL	Cylinder
AXIS	Axis
L	Left
R	Right

Date of preparation:

2015-03

*a trademark of Novartis